




FlexArc[™] MicroCatheter

FDA 510(K) No. K243534

CE MDR 845697

NMPA 20223030593

DEVELOPING INNOVATIVE
TECHNOLOGIES TO OPTIMIZE
VASCULAR TREATMENT




● 80 cm hydrophilic coating.


● Balanced shaft transitions and robust column strength minimize tip movement, providing reliable retrieval of treatment devices.

● Stenghtened stainless steel braiding leads to improved maneuverability and pushability.

● 1 mm tip ensures safe operation inside blood vessels while better selecting blood vessels.



Peel- Able Sheath



Shape Needles

Product Code	ID [inch/mm]	Distal OD [mm]	Proximal OD [mm]	Working length [cm]
200017150*	0.017" /0.43	0.66	0.81	150
200017153	0.017" /0.43	0.66	0.81	153
200017155*	0.017" /0.43	0.66	0.81	155
200021150*	0.021" /0.53	0.81	0.91	150
200021153	0.021" /0.53	0.81	0.91	153
200021155*	0.021" /0.53	0.81	0.91	155
200027150*	0.027" /0.69	0.93	0.93	150
200027155*	0.027" /0.69	0.93	0.93	155

* FDA 510(K) clearance



Zenith Vascular Scitech

TEL: 86-0512-69138730

E-mail: International@zenithmedtech.com

Website: Zenithmedtech.com