


Zenith Delivery Catheter

NMPA 20243031800



Developing Innovation Technologies To Optimize Vascular Treatment

● New Choice for Radial Access

- 6F&7F radial-specific delivery catheters
- Optimized segment design for radial access, offering stronger support in the aortic arch.

● Better Choice for Femoral Access

- 8F delivery catheter with the largest inner lumen (0.090 inch) on the market.
- Easily compatible with multiple microcatheter systems, and allows smooth device advancement with some extra lumen space.

● Compatible with Innovative Procedures

- A wide range of sizes: OD from 6 to 10F and length from 70 to 105 cm.
- Provides comprehensive access support for innovative procedures.



OD	Working Length				
	70cm	80cm	90cm	100cm	105cm
6F	213072070	213072080	213072090	213072100	213072105
	213072070MP	213072080MP	213072090MP	213072100MP	213072105MP
7F	213082070	213082080	213082090	213082100	213082105
	213082070MP	213082080MP	213082090MP	213082100MP	213082105MP
8F	213090070	213090080	213090090	213090100	213090105
	213090070MP	213090080MP	213090090MP	213090100MP	213090105MP
9F	213100070	213100080	213100090	213100100	213100105
	213100070MP	213100080MP	213100090MP	213100100MP	213100105MP
10F	/	213111080	213111090	213111100	213111105
	/	213111080MP	213111090MP	213111100MP	213111105MP



Zenith Vascular Scitech

TEL: 86-0512-6913873

E-mail: Information@zenithmedtech.com

Website: Zenithmedtech.com